



Tel: 1-909-868-1800  
Fax: 1-909-868-1898  
www.zenithlabo.com  
Email: info@zenithlabo.com  
1979 W Holt Ave, Pomona, CA 91768

# MONOCULAR BIOLOGICAL MICROSCOPE



MODEL	XSP-1CA
Total Magnification	40X~1600X
Eyepieces	WF10X, WF16X
Objectives	4X, 10X, 40Xs, 100Xs
Eyepieces Tube	Binocular 45° inclined Tube Rotatable 360°
Illumination	<ul style="list-style-type: none"><li>● 6V/20W Halogen Lamp</li><li>● 20W Incandescent</li></ul>
Coarse Adjustment Range	14mm
Working Stage	<ul style="list-style-type: none"><li>● Double Layer Mechanical Stage:110x125mm</li><li>● Moving Range 60mmx20mm</li></ul>
Condenser	NA1.25 Abee Condenser With Changeable Diaphragm, Moving Scale And Filter
Size	36x26x21cm/14.2x10.3x8.3in
Gross Weight	3.5kg/8lb
Net Weight	3kg/7lb