



Tel: 1-909-868-1800

Fax: 1-909-868-1898

www.zenithlabo.com

Email: info@zenithlabo.com

1979 W Holt Ave, Pomona, CA 91768

HIGH SPEED CENTRIFUGE HC - 16L

- Strong plastic case with metal lid.
- Brushless motor, setting by speed, RCF and time.
- LCD display speed, RCF and time.
- Over speed and imbalance alarm.
- Touch panel control, automatic stop when lid open.
- Error code display, speed up and down adjustment.
- Suitable for researches on genes and protein nucleic acid.



MODEL	HC- 16L
Max RPM	16000 RPM. Stepless speed adjustment, LED display speed, RCF and time.
MAX RCF	17000xg
TIME	0-99 min
Rotor Capacity	1.5/2mlx12 (Standard)
Voltage	110V 60Hz/ 220V 50Hz
Dimensions (mm)	320x340x280
Weight	9 kg

Model	Ordering number	Max speed	Rotor capacity	Max RCF
HC-16L	A160212	16000	1.5/2mlx12	17000xg
	A160506	16000	5mlx6	17000xg
	A160012	16000	0.5mlx12	17000xg
	A130508	13000	5mlx8	11400xg
	A120224	12000	1.5/2mlx24	15100xg